

A COMPLETE SAFETY SYSTEM.

FLASH thermosealer



FLASH, thermosealer

- Ergonomic design
- Wall mountable
- Compact size
- Practical and easy to use
- Quick sealing
- Packing support (optional)



FLASH

	U.M.	
Packet tray		optional
Overall dimensions W x H x D	mm	445 x 220 x 330
Electrical voltage and frequency	V - Hz	230 - 50 (single - phase AC)
Absorbed maximum power	W	100
Seal width	mm	10
Seal length	mm	250
Seal duration	s	3
Weight (mass)	kg	5,2

Technical data are subject to change without notice