

SONICLEAN



SONICLEAN

- Available in two sizes, 3 and 6 liters, completely in stainless steel
- Lid and basket included
- Programmable timer from 0 to 30 minutes, with the possibility of continuous operation
- Programmable heater from 0°C to 80°C
- SWEEP function for the perfect ultrasound distribution in the cleaning liquid
- BOOST function to increase the ultrasound power by approx. 25 %
- DEGAS function for the efficient degassing of the cleaning liquid
- Quick drain system on the back, standard on SoniClean 6 model
- Durable and reliable
- Easy to clean

SONICLEAN

	M.U.	SoniClean 3	SoniClean 6
Tank capacity	l	2,75	5,75
Drain system			•
Lid		•	•
Basket		•	•
Becker holder		optional	optional
Becker		optional	optional
Tank internal dimensions W x H x D	mm	240 x 137 x 100	300 x 151 x 150
Basket internal dimensions W x H x D	mm	198 x 106 x 50	255 x 115 x 75
Unit external dimensions W x H x D	mm	300 x 179 x 214	365 x 186 x 264
Electrical voltage and frequency	V - Hz	230 - 50/60 (single-phase AC)	
Absorbed maximum power	W	280	550
Weight (mass)	kg	3,3	5,1
Ultrasonic frequency	kHz	37	

• Standard equipment

Technical data are subject to change without notice